

EVAC+CHAIR®

3000 Marcus Ave, Suite #3E6, Lake Success, NY 11042-1012

Tel: (516) 502-4240 Fax: (516) 327-8220

www.evac-chair.com sales@evac-chair.com

How would you get these people down the stairs, if the Lifts or Elevators do not work due to power failure or even worse, a fire?

The **ResQme*** is the easy,
fast and safe way to
Evacuate



1. Transfer the person on to the ResQme



2. Slide down stairs



3. Controlled descent

“Slide to Safety”



Nursing Homes




Apartments



Hotels

- Pulling Straps at both ends
- Head pillow for comfort and protection
- Will slide on Carpet, Tiles, Vinyl and Concrete
- The ResQme can be pulled by one or up to four people
- Can be used for children, small adults and people up to 6' 6" tall
- Strong 4" wide Velcro Straps to secure the person being rescued
- Wrist/Arm strap to stop hands and arms from flaying out of control
- Base board to protect the person from hitting door frames etc., when going around corners

Weight Tested to: 660 lbs.
S.W.L. of 440lbs.
In accordance with
AS ISO 10535 – 2002



Hospitals



Shopping Centres



Skyscrapers